


# Technical Specification

## O2-G2 Oxygen Sensor Miniature Size



Figure 1 O2-G2 Schematic Diagram



Top View

Bottom View

Side View

<b>PERFORMANCE</b>	Output Response time Zero current	$\mu\text{A}$ @ 22°C, 20.9% O <sub>2</sub> t <sub>90</sub> (s) from 20.9% to 0% O <sub>2</sub> (47Ω) $\mu\text{A}$ @ 99.99% N <sub>2</sub> , 22°C	30 to 42 < 18 < 2
<b>LIFETIME</b>	Output drift Operating life	% change in output @ 3 months months until 85% original output in 20.9% O <sub>2</sub>	< 2 > 24

### ENVIRONMENTAL

Humidity Sensitivity	% O <sub>2</sub> change: 0% to 95% rh @ 40°C	< 0.7
Pressure sensitivity	(% change of output)/(% change of pressure) @ 20kPa	< 0.1
CO <sub>2</sub> sensitivity	% change in output / % CO <sub>2</sub> @ 5% CO <sub>2</sub>	< 0.1
Output at -20°C	%output/output at 20°C in 20.9% O <sub>2</sub>	87 to 93
Output at +50°C	%output/output at 20°C in 20.9% O <sub>2</sub>	103 to 107

### KEY SPECIFICATIONS


Temperature range	°C	-30 to 55
Pressure range	kPa	80 to 120
Humidity range	% rh continuous (0 to 99% rh short term)	5 to 95
Storage period	months @ 3 to 20°C (store in sealed pot)	6
Load resistor	Ω (recommended)	47 to 100
Weight	g	< 7

**NOTE:** all sensors are tested at ambient environmental conditions, with 47 ohm load resistor, unless otherwise stated. As applications of use are outside our control, the information provided is given without legal responsibility. Customers should test under their own conditions, to ensure that the sensors are suitable for their own requirements.

# Technical Specification

## O2-G2 Performance Data


Figure 2 Temperature Dependence in Air



This graph shows the variation in sensitivity caused by changes in temperature.

All capillary oxygen sensors will show some variation in signal output with temperature and the typical response of an O2-G2 is shown.


Figure 3 Pressure Step Performance



Step changes in pressure can cause a temporary signal transient. Positive pressure gives a output signal increase whilst negative pressure causes the output signal to decrease.

Typical transient response for an O2-G2 sensor exposed to a 10kPa pressure pulse is shown.

Figure 4 Linearity



Mass flow oxygen sensors generate a non-linear current with increasing oxygen concentration:  
current =  $k * \log(1/(1-C))$ .

When plotted on a linear graph, figure 4 shows that the non-linearity is very repeatable and can be corrected in software to the required accuracy.